


Hydrangea as Cut Flower

Prepared by Maher Jabado

- 1) You need to have:
- -Fresh Hydrangea (better cut in the morning after good watering one day earlier)
- -500 ml warm water
- -Liquid dish detergent
- -Sugar
- -White Vinegar
- -Mouth wash antiseptic (i.e Listerine)
- -Hammer
- 2) -Make a fresh cut (diagonal) of the Hydrangea stem; I prefer to have the cut made while stem is soaked in water
- 3)-SOAK the stem completely in COLD WATER and keep there for 2 hours
- 4) -Prepare the following preserved solution:
- a) Inside the 500 ml warm solution place 5 drops of







- (i.e antiseptic Listerine)
- 5) Shake well the above mix
- 6) Remove the lower Hydrangea leaves that could come under water level in the vase
- 7) Use the hammer and hit a few cm at the bottom of the stem

